

01

Diagnotics-Anesthesia

Diagnostik- Anästhesie
Diagnostico - Anestesia
Diagnostic- Anesthésie
Diagnostica - Anestesia



Babinsky
CFS-01-100-22
 22 cm - 8 3/4"

Rabiner
CFS-01-101-23
 23 cm - 9 1/4"

Queen's Hammer
CFS-01-108-35
 35 cm - 13 3/4"
 Plastic handle



Wartenberg
CFS-01-102-19
19 cm - 7 1/2"



Buck
CFS-01-103-18
18 cm - 7"



Traube
CFS-01-104-18
18 cm - 6 1/4"



Taylor
CFS-01-105-18
18 cm - 7"
CFS-01-105-20
20 cm - 8"



Berliner
CFS-01-106-20
20 cm - 8"



Dejerine
CFS-01-107-20
20 cm - 8"

Diagnostic Sets

Trousses a diagnose
Juegos para diagnose
Astucci diagnostici
Diagnostik-Bestecke

Set - 01

CFS-01-108-00

Ophthalmoscope set
Ophthalmoskop-Besteck
Ophthalmoscopio
Trousse à ophtalmoscopie
Ofthalmoscopio



Set - 02

CFS-01-109-00

Universal set
Universal-Besteck
Estuche universal
Trousse universelle
Astuccio universale



Diagnostic Sets

Trousses a diagnose
Juegos para diagnose
Astucci diagnostici
Diagnostik-Bestecke



Set-03

CFS-01-110-00

Combination set
Kombinations-Besteck
Estuche combinado
Trousse combinée
Astuccio combinato



Set-04

CFS-01-111-00

Otoscope set
Otoskop-Besteck
Otoscopio
Trousse à otoscopie
Otoscopio





Pinard
CFS-01-112-18
 18 cm - 7"
 Wood
 Holz
 Madera
 Bois
 Legno




Pinard
CFS-01-113-15
 15 cm - 6"
 Aluminium
 Aluminium
 Aluminio
 Aluminium
 Alluminio



Duplex
CFS-01-114-04 4.7cm - 1 3/4"
 for adults
 für Erwachsene
 para adultos
 pour adultes
 per adulti



Baby-Duplex
CFS-01-115-15 4.7cm - 1 3/4"
 for children
 für Kinder
 para niños
 pour enfants
 per bambini



Anestophon
CFS-01-116-04 4.7cm - 1 3/4"
 for nurses, simple
 für Schwestern, einfach
 para enfermeras, simple
 pour infirmières, simple
 per infermiere, semplice





Ziegler

CFS-01-117-09 9 cm - 3 1/2" Ø

Head mirror with black plastic head band
 Stirnspiegel mit schwarzem Plastikband
 Espejo frontal con banda de plastica negra
 Miroir frontal avec ruban en plastique noir
 Specchio frontale con fascia di plastica nera

Ziegler

CFS-01-118-09 9 cm - 3 1/2" Ø

Head mirror only, with rubber coating
 Reflexspiegel allein, mit Gummifassung
 Espejo frontal solo, con guarnición de goma
 Miroir frontal seul, avec monture en caoutchouc
 Solo specchio frontale, con guarnizione di gomma

CFS-01-119-09 9 cm - 3 1/2" Ø

Head mirror with white plastic head band
 Stirnspiegel mit weissem Plastikband
 Espejo frontal con banda de plastica blanca
 Miroir frontal avec ruban en plastique blanc
 Specchio frontale con fascia di plastica bianca

CFS-01-120-00

Black plastic head band
 Plastikband, schwarz
 Banda de plastica, negra
 Ruban en plastique noir
 Fascia di plastica, nera

CFS-01-121-00

White plastic head band
 Plastikband, weiss
 Banda de plastica, blanca
 Ruban en plastique blanc
 Fascia di plastica, bianca

Cold-light head lamp
 Kaltlicht-Stirnlampe
 Lámpara frontal de luz fria
 Eclairer frontal de lumière froide
 Lampada frontale a luce fredda



Binner

CFS-01-122-00

Head light, cold-light, adjustable focus,
 working distance 300 mm, field of illumination 70 –120 mm Ø

Kaltlicht-Stirnlampe mit verstellbarer Fokuserlinse,
 Arbeitsabstand 300mm, Leuchtfeldgröße 70 – 120 mm Ø

Lámpara frontal, luz fria, focus ajustable,
 distancia de trabajo 300mm, campo de iluminación 70 – 120mm Ø

Eclairer frontal, lumière froide, focus ajustable,
 distance de travail 300 mm, champ d'illumination 70 – 120mm Ø

Lampada frontale, luce fredda, con lente a focale variabile,
 distanza di lavoro 300 mm, campo d'illuminazione 70 – 120mm Ø





Clar
CFS-01-123-00 10 cm Ø, 6 V

Head light complete
Stirnlampe komplett
Lámpara frontal, completa
Eclaireur frontal complet
Lampada frontale, completa




CFS-01-124-00 10 cm Ø, 6 V

Head light with cable, frame and mirror, without head band
Stirnlampe mit Kabel, Fassung und Spiegel, ohne Stirnreifen
Lámpara frontal, con cable, armadura y espejo, sin banda
Eclaireur frontal avec cable, douille et miroir, sans ruban
Lampada frontale con cavo, supporto e specchio, senza fascia




Mini-Clar
CFS-01-125-00 5 cm Ø, 6 V

complete
komplett
completa
complet
completa



Clar
CFS-01-126-05 5 cm - 2" Ø




CFS-01-126-10 10 cm - 2" Ø

Spiegel aus Glas, ohne Fassung
Head mirror glass, without frame
Miroir frontal en vitre, sans douille
Espejo frontal de cristal, sin armadura
Vetro specchio, senza supporto



Clar
CFS-01-127-00

Spezial-Gelenkverbindung
Special socket joint
Assemblage articulé spécial
Conexión
Giunto di connessione speciale



Clar
CFS-01-128-00

Stirnreifen allein
Head band only
Ruban seul
Banda sola
Solo fascia frontale



CFS-01-129-00 6 V

Glühbirne
Bulb
Ampoule
Bombilla
Lampadina



Transformer – stepless adjustable
Transformator – stufenlos einstellbar
Transformador – regulación continua
Transformateur – à réglage continuelle
Trasformatore – regolabile in continuità

CFS-01-130-00

Input 220 V/ 50 Hz~
Output 0–10 V/ 3 W=

Binocular

CFS-01-131-00

2.5 times enlargement, 340 mm working distance, field of vision 90 mm Ø, complete.
Pupillary distance adjustable from 48 mm–77 mm. Individual adjustable angle of inclination.
Weight of the optic 42g.

2.5 fache Vergrößerung, 340 mm Arbeitsabstand, Sichtfeld > 90 mm Ø, komplett.
Pupillenabstand einstellbar von 48 mm–77 mm. Neigungswinkel individuell justierbar.
Gewicht der Optik 42g.

2.5 veces aumento, distancia de trabajo 340 mm, campo de visión >90mm Ø, completo.
Distancia interpupilar ajustable de 48 mm–77 mm. Angulo de inclinación ajustable individual.
Peso de la optica 42g.

2.5 fois grossissement, distance de travail 340 mm, champ de vision >90mm Ø, complet.
Ecart interpupillaire ajustable de 48 mm–77 mm. Angle d'inclinaison ajustable individuellement.
Poids de l'optique 42g.

Ingrandimento: 2.5 volte, distanza di lavoro 340 mm, campo di visione >90mm Ø, completo.
Distanza interpupillare regolabile da 48 mm–77 mm. Angolo di inclinazione regolabile individualmente.
Peso dell'ottica 42g.



CFS-01-121-00

CFS-01-132-00

- Kopfband
- Head band
- Ruban
- Banda
- Fascia

- Sphygmomanometers
- Blutdruckmessgeräte
- Esfigmomanómetros
- Tensiomètres
- Sfigmomanometri

**Manometer
Fig. 1**

CFS-01-133-00

- Calibrated
- Geeicht
- Calibrado
- Calibré
- Calibrato



**Manometer
Fig. 2**

CFS-01-134-00


- non calibrated
- ungeeicht
- no calibrado
- non calibré
- non calibrato



**Manometer
Fig. 3**

CFS-01-135-00

- Mercury blood pressure manometer, non calibrated
- Quecksilber-Blutdruckmessapparat, ungeeicht
- Manómetro de mercurio, no calibrado
- Tensiomètre à mercure, non calibré
- Manometro a mercurio, non calibrato





Collin
 Fig. 1

CFS-01-136-12

12.5 cm - 5"
 for children, kg-scale
 für Kinder, kg-Einteilung
 para niños, graduación en kg
 pour enfants, échelle en kg
 per bambini, graduazione in kg

Fig. 2

CFS-01-137-12

12.5 cm - 5"
 for adults, kg-scale
 für Erwachsene, kg-Einteilung
 para adultos, graduación en kg
 pour adultes, échelle en kg
 per adulti, graduazione in kg

Fig. 3

CFS-01-138-12

12.5 cm - 5"
 for adults, lbs-scale
 für Erwachsene, lbs-Einteilung
 para adultos, graduación en lbs
 pour adultes, échelle en lbs
 per adulti, graduazione in lbs



Bertillon

CFS-01-139-32

32 cm - 12 1/2"



Moeltgen

CFS-01-140-24

24 cm 9 1/2"





Magda
 Fig. 1

CFS-01-141-10	10 cm - 4"
CFS-01-141-15	15 cm - 6"
CFS-01-141-20	20 cm - 8"
CFS-01-141-30	30 cm - 12"
CFS-01-141-50	50 cm - 20"



Magda
 Fig. 2

CFS-01-142-50	50 cm - 20"
X-ray ruler	
Röntgen-Messstab	
Regla para rayos X	
Règle rayons X	
Misuratore per raggi X	



Vernier

CFS-01-143-21	CFS-01-144-21	21 cm - 8 1/4"
Aluminium	Steel	
Aluminium	Stahl	
Aluminio	Acero	
Aluminium	Acier	
Alluminio	Acciaio	



Castroviejo

CFS-01-145-08	8 cm - 3 1/4"
---------------	---------------



Gauge

CFS-01-146-15	CFS-01-147-15	150 cm - 60"
Steel	Wax cloth	
Stahl	Wachs	
Acero	Cerada	
Acier	Toile cirée	
Acciaio	Cera	





<p>Welzer CFS-01-152-16 16 cm - 6 ¼"</p> <p>Wood Holz Madera Bois Legno</p>	<p>Wavy CFS-01-153-14 14.5 cm - 5 ¾"</p>	<p>Kroenig CFS-01-154-15 15 cm - 6"</p>	<p>Mayo CFS-01-155-17 17 cm - 6 ¾"</p>
<p>Bruenings CFS-01-156-16 19 cm - 7 ½"</p>	<p>Frenzel CFS-01-157-16 16.5 cm - 6 ½"</p>	<p>Lack CFS-01-130-13 13mm CFS-01-130-19 19mm CFS-01-130-22 22mm</p>	<p>Wieder CFS-01-158-15 15 cm - 6"</p> <p>CFS-01-158-14 14 cm - 5 ½"</p>



Bosworth

CFS-01-159-14

14.5 cm - 5 3/4"



Tobold

CFS-01-160-14

14.5 cm - 5 3/4"



Hartmann

CFS-01-161-15

15 cm - 6"



Andrew

CFS-01-162-11

11 cm - 4 3/4"

Tongue holding forceps
 Zungenzangen
 Pinzas tiralingua
 Pincas tire langue
 Pinze tirilingua



Young

CFS-01-163-17

17 cm - 6 3/4"

CFS-01-163-00

- 1 Pair rubber-inserts
- 1 Paar Ersatz Gummi-Einsätze
- 1 Juego de insertos de goma
- 1 Paire des inserts en caoutchouc
- 1 Paio d'inserti di gomma



Collin

CFS-01-164-16 16 cm - 6 1/4"

CFS-01-164-17 17 cm - 6 3/4"

CFS-01-164-19 19 cm - 7 1/2"



Cheek Retractors

Wangenhalter

Separador de mejillas

Ecarteur des joues

Divaricatori per guancia



Sternberg

CFS-01-165-12 12 cm - 4 3/4"

CFS-01-165-14 14 cm - 5 1/2"

CFS-01-165-16 16 cm - 6 1/4"

13 Catheter forceps
Katheterzangen
Pinzas para catéteres
12 Pincas pour cathéters
Pinze per catetere



Magill

CFS-01-166-17

17 cm - 6"
Infant



CFS-01-166-20

20 cm - 8"
Child



CFS-01-166-25

25 cm - 9"
Adult



! The English Macintosh is similar to the conventional Macintosh except that the flange is curved and lower at the handle end. Two studies have proved that this blade provided better results than other commonly used blades.



Fig. 00



Fig. 0



Fig. 1



Fig. 2



Fig. 3



Fig. 4




Fig. 5

Code	Description	Fig.	Size
CFS-01-281-00	Premature	00	68 mm
CFS-01-281-0	Neonate	0	80 mm
CFS-01-281-1	Infant	1	92 mm
CFS-01-281-2	Child	2	105 mm
CFS-01-281-3	Medium Adult	3	135 mm
CFS-01-281-4	Large Adult	4	155 mm
CFS-01-281-5	X-Large	5	175 mm





! The Miller is one of the most popular blades. The tongue is straight with a slight upward curve near the tip. In cross-section the flange, web, and tongue form a C with the top fattened. Placing the bulb on the right side allows it to be protected by the flange. Several modifications have been described. In comparing the Macintosh and Miller blades, it was found that force, head extension and cervical spine movement were less with the Miller.



Code	Description	Fig.	Size
CFS-01-282-00	Premature	00	65 mm
CFS-01-282-0	Neonate	0	78 mm
CFS-01-282-1	Infant	1	102 mm
CFS-01-282-2	Child	2	154 mm
CFS-01-282-3	Medium Adult	3	195 mm
CFS-01-282-4	Large Adult	4	205 mm

- ! Specially designed and constructed with high quality stainless steel AISI 304.
- High quality German optical fibers loaded with 4.5mm bundle for best optimum light transmission.
- Exceptionally high light intensity with LED illumination (8000-10000 LUX).
- With this integrated FO light conductor, there are no extra cavities or light path to become contaminated.
- Smooth and low profile blade design which gives no chance of damaging teeth and upper part of mouth.
- Integrated light guide blade design for quick cleaning and sterilization without extra effort for connecting screw.
- Compatible with all handles of ISO 7376 (Green System).
- Autoclaving at 134 C /5 min. maximum light output with 3000-4000 cycles.



Fig. 0



Fig. 1



Fig. 2



Fig. 3



Fig. 4



Fig. 5

Code	Description	Fig.	Size
CFS-01-283-0	Neonate	0	80 mm
CFS-01-283-1	Infant	1	93 mm
CFS-01-283-2	Child	2	114 mm
CFS-01-283-3	Medium Adult	3	135 mm
CFS-01-283-4	Large Adult	4	155 mm
CFS-01-283-5	X-Large	5	175 mm



- ! Specially designed and constructed with high quality stainless steel AISI 304.
- High quality German optical fibers loaded with 4.5mm bundle for best optimum light transmission.
- Exceptionally high light intensity with LED illumination (8000-10000 LUX).
- With this integrated FO light conductor, there are no extra cavities or light path to become contaminated.
- Smooth and low profile blade design which gives no chance of damaging teeth and upper part of mouth.
- Integrated light guide blade design for quick cleaning and sterilization without extra effort for connecting screw.
- Compatible with all handles of ISO 7376 (Green System).
- Autoclaving at 134 C /5 min. maximum light output with 3000-4000 cycles.



Code	Description	Fig.	Size
CFS-01-284-00	Premature	00	65 mm
CFS-01-284-0	Neonate	0	78 mm
CFS-01-284-1	Infant	1	102 mm
CFS-01-284-2	Child	2	154 mm
CFS-01-284-3	Medium Adult	3	195 mm
CFS-01-284-4	Large Adult	4	205 mm



! These blades have a flexible tip that is controlled by a lever attached to the proximal end of the blade. When the lever is pushed toward the handle, the tip of the blade is flexed. This laryngoscope may be helpful when a difficult intubation is encountered. It may be especially useful in patients with minimal neck movement. External laryngeal pressure may further improve the view. Even if the blade does not improve the view of the larynx, it may improve the likelihood of successful intubation by elevating the epiglottis.




Fig. 1



Fig. 2



Fig. 3

Code	Description	Fig.	Size
CFS-01-285-2	Child	2	105 mm
CFS-01-285-3	Medium Adult	3	130 mm
CFS-01-285-4	Large Adult	4	155 mm





Fig. 2



Fig. 3



Fig. 4

! These blades have a flexible tip that is controlled by a lever attached to the proximal end of the blade. When the lever is pushed toward the handle, the tip of the blade is flexed. This laryngoscope may be helpful when a difficult intubation is encountered. It may be especially useful in patients with minimal neck movement. External laryngeal pressure may further improve the view. However, it will not always improve the view and may worsen it. Even if the blade does not improve the view of the larynx, it may improve the likelihood of successful intubation by elevating the epiglottis. In some patients in whom the view is satisfactory with the Macintosh blade, use of the Flexible tip blade may make the view worse.

Code	Description	Fig.	Size
CFS-01-286-2	Child	2	108 mm
CFS-01-286-3	Medium Adult	3	130 mm
CFS-01-286-4	Large Adult	4	155 mm

Battery Handles for Fiber Optic Laryngoscope

Batteria Griffe fur Fiber OPTic Laryngoskop

Manijas de la bateria para el sistema verde laringoscopio de fibre optique

Poingnees de batterie pour le systeme vert laryngoscope a fibre optique

Batteria maniglie per sistema di fibre ottiche laringoscopio Verde



- ! Diamond knurled barrel gives strong grip to the user.
- LED brighter illumination with 60,000 LUX.
- LED operating time over 50,000 hrs.
- 2 hrs operating time at maximum illumination.
- Viewing angle at 130 with color temperature color temperature 6500K.
- Conforms to all Green System blades as per ISO 7376 standard.
- Simple change of batteries at the base of the handle.
- Crafted from high quality chrome plated brass.
- Battery usage is minimized due to the low power consumption.
- Fully autoclavable @ 134 C for 5 minutes.



Battery Handles Fiber Optic Xenon, Stainless Steel

Code	Description	Battery	Size
CFS-01-287-01	Infant	1xLithium 123 Size	122 mm
CFS-01-287-02	Stubby	2xAA Size	160 mm
CFS-01-287-03	Pediatric	2xAA Size	125 mm
CFS-01-287-04	Medium	2xC Size	160 mm

Battery Handles Fiber Optic Xenon, Brass Chrome Plated

Code	Description	Battery	Size
CFS-01-288-01	Infant	1xLithium 123 Size	122 mm
CFS-01-288-02	Stubby	2xAA Size	160 mm
CFS-01-288-03	Pediatric	2xAA Size	125 mm
CFS-01-288-04	Medium	2xC Size	160 mm

Battery Handles Fiber Optic LED, Stainless Steel

Code	Description	Battery	Size
CFS-01-289-01	Infant	1xLithium 123 Size	122 mm
CFS-01-289-02	Stubby	2xAA Size	160 mm
CFS-01-289-03	Pediatric	2xAA Size	125 mm
CFS-01-289-04	Medium	2xC Size	160 mm

Battery Handles Fiber Optic LED, Brass Chrome Plated

Code	Description	Battery	Size
CFS-01-290-01	Infant	1xLithium 123 Size	122 mm
CFS-01-290-02	Stubby	2xAA Size	160 mm
CFS-01-290-03	Pediatric	2xAA Size	125 mm
CFS-01-290-04	Medium	2xC Size	160 mm






Battery Handles Fiber Optic LED, Stainless Steel

Code	Description	Battery	Size
CFS-01-291-01	Infant	1xLithium 123 Size	122 mm
CFS-01-291-02	Stubby	2xAA Size	160 mm
CFS-01-291-03	Pediatric	2xAA Size	125 mm
CFS-01-291-04	Medium	2xC Size	160 mm



Replacement Lamps

Code	Type	Voltage
CFS-01-292-01	Xenon	2.5 V
CFS-01-292-02	LED	3.0 V
CFS-01-292-03	Halogen	2.5 V
CFS-01-292-04	Halogen	3.5 V



Replacement Fiber Optic Light Guides, English Profile (4 mm)

Code	Description	For Blades
CFS-01-296-00	Miller Blade	00
CFS-01-296-0	Miller Blade	0
CFS-01-296-1	Miller Blade	1
CFS-01-296-2	Miller Blade	2
CFS-01-296-3	Miller Blade	3
CFS-01-296-4	Miller Blade	4
CFS-01-296-00	Macintosh Blade	00
CFS-01-296-0	Macintosh Blade	0
CFS-01-296-1	Macintosh Blade	1
CFS-01-296-2	Macintosh Blade	2
CFS-01-296-3	Macintosh Blade	3
CFS-01-296-4	Macintosh Blade	4



Macintosh Fiber Optic Xenon Laryngoscope Sets, English Profile

Code	Description
CFS-01-298-01	Set of 4 Blades 1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
CFS-01-298-03	Set of 5 Blades 0,1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
CFS-01-298-04	Set of 6 Blades 0,1,2,3,4,5 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.



Miller Fiber Optic Xenon Laryngoscope Sets, English Profile

Code	Description
CFS-01-299-01	Set of 4 Blades 1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
CFS-01-299-03	Set of 5 Blades 0,1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
CFS-01-299-04	Set of 6 Blades 0,1,2,3,4,5 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.



Flexible Fiber Optic Xenon Laryngoscope Sets, English Profile

Code	Description
CFS-01-300-01	Set of 2 Blades 3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.

Flexible pro Fiber Optic Xenon Laryngoscope Sets, German Profile


Code	Description
CFS-01-301-01	Set of 2 Blades 3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.



Macintosh Fiber Optic LED Laryngoscope Sets, English Profile

Code	Description
CFS-01-302-01	Mac FO LED Laryngoscope Set of 4 Blades 1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
CFS-01-302-02	Mac FO LED Laryngoscope Set of 5 Blades 0,1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
CFS-01-302-03	Mac FO LED Laryngoscope Set of 6 Blades 0,1,2,3,4,5 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.





Miller Fiber Optic LED Laryngoscope Sets, English Profile

Code	Description
CFS-01-303-01	Miller FO LED Laryngoscope Set of 4 Blades 1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
CFS-01-303-02	Miller FO LED Laryngoscope Set of 5 Blades 0,1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
CFS-01-303-03	Miller FO LED Laryngoscope Set of 6 Blades 00,0,1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.



Flexible Fiber Optic Led Laryngoscope Sets, English Profile

Code	Description
CFS-01-304-01	MacFlex FO LED Laryngoscope Set of 2 Blades 3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.

Flexible pro LED Laryngoscope Sets, German Profile

Code	Description
CFS-01-305-01	Set of 2 Blades 3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.

I The English Macintosh is similar to the conventional Macintosh except that the flange is curved and lower at the handle end. Two studies have proved that this blade provided better results than other commonly used blades.



Fig. 00



Fig. 0



Fig. 1



Fig. 2



Fig. 3



Fig. 3.5



Fig. 4



Fig. 5

Code	Description	Fig.	Size
CFS-01-306-00	Premature	00	68 mm
CFS-01-306-0	Neonate	0	80 mm
CFS-01-306-1	Infant	1	95 mm
CFS-01-306-2	Child	2	105 mm
CFS-01-306-3	Medium Adult	3	135 mm
CFS-01-306-35	Ext. Medium Adult	3.5	145 mm
CFS-01-306-4	Large Adult	4	155 mm
CFS-01-306-5	X-Large	5	175 mm



! The Miller is one of the most popular blades. The tongue is straight with a slight upward curve near the tip. In cross-section the flange, web, and tongue form a C with the top fattened. Placing the bulb on the right side allows it to be protected by the flange. Several modifications have been described. In comparing the Macintosh and Miller blades, it was found that force head extension and cervical spine movement were less with the Miller.



Fig. 00



Fig. 2



Fig. 0



Fig. 3



Fig. 1



Fig. 4

Code	Description	Fig.	Size
CFS-01-307-00	Premature	00	65 mm
CFS-01-307-0	Neonate	0	78 mm
CFS-01-307-1	Infant	1	102 mm
CFS-01-307-2	Child	2	154 mm
CFS-01-307-3	Medium Adult	3	195 mm
CFS-01-307-4	Large Adult	4	205 mm

! These blades have a flexible tip that is controlled by a lever attached to the proximal end of the blade. When the lever is pushed toward the handle, the tip of the blade is flexed. This laryngoscope may be helpful when a difficult intubation is encountered. It may be especially useful in patients with minimal neck movement. External laryngeal pressure may further improve the view. Even if the blade does not improve the view of the larynx, it may improve the likelihood of successful intubation by elevating the epiglottis.



Code	Description	Fig.1	Size
CFS-01-308-2	Child	2	105 mm
CFS-01-308-3	Medium Adult	3	130 mm
CFS-01-308-4	Large Adult	4	155 mm



Battery Handles for Fiber Optic Laryngoscope

Batteria Griffe fur Fiber OPTic Laryngoskop

Manijas de la bateria para el sistema verde laringoscopio de fibre optique

Poingnees de batterie pour le systeme vert laryngoscope a fibre optique

Batteria maniglie per sistema di fibre ottiche laringoscopio Verde




- ! Diamond knurled barrel gives strong grip to the user.
- Conforms to all conventional blades as per ISO 7376 standard.
- Simple change of batteries at the base of the handle.
- Crafted from high quality chrome plated brass.
- Economical solution for intubation as compare to green system.
- Fully autoclavable @ 134 C for 5 minutes.

Battery Handles, Stainless Steel

Code	Description	Battery	Size
CFS-01-309-01	Infant	1xLithium 123 Size	110 mm
CFS-01-309-02	Stubby	2xAA Size	115 mm
CFS-01-309-03	Pediatric	2xAA Size	145 mm
CFS-01-309-04	Medium	2xC Size	150 mm

Battery Handles, Brass Chrome Plated

Code	Description	Battery	Size
CFS-01-310-01	Infant	1xLithium 123 Size	110 mm
CFS-01-310-02	Stubby	2xAA Size	115 mm
CFS-01-310-03	Pediatric	2xAA Size	145 mm
CFS-01-310-04	Medium	2xC Size	150 mm



Replacement Lamps

Code	Description	Type	Voltage
CFS-01-311-00	Small	Xenon	2.5 V
CFS-01-311-00	Adult	LED	3.0 V
CFS-01-311-00	Small	Halogen	2.5 V
CFS-01-311-00	Adult	Halogen	3.5 V



Macintosh Conventional Laryngoscope Sets, English Profile

Code	Description
CFS-01-315-01	Mac Conventional Laryngoscope Set of 4 Blades 1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
CFS-01-315-02	Mac Conventional Laryngoscope Set of 5 Blades 0,1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
CFS-01-315-03	Mac Conventional Laryngoscope Set of 6 Blades 00,0,1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.



Miller Conventional Laryngoscope Sets, English Profile

Code	Description
CFS-01-316-01	Miller Conventional Laryngoscope Sets Set of 4 Blades 1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp in Plastic Case.
CFS-01-316-02	Miller Conventional Laryngoscope Sets Set of 5 Blades 0,1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.
CFS-01-316-03	Miller Conventional Laryngoscope Sets Set of 6 Blades 00,0,1,2,3,4 with 1 Medium Handle, 1 Replacement Lamp, in Plastic Case.



Inka
CFS-01-192-05 x = 17.2 cm

Fig. 5
 Stainless steel
 Rostfreier Edelstahl
 Acero inoxidable
 Acier inoxydable
 Acciaio inossidabile



Inka
CFS-01-193-06 x = 18.2 cm

Fig. 6
 Stainless steel
 Rostfreier Edelstahl
 Acero inoxidable
 Acier inoxydable
 Acciaio inossidabile



Inka
CFS-01-194-07 x = 22.2 cm

Fig. 7
 Rostfreier Edelstahl
 Stainless steel
 Acier inoxydable
 Acero inoxidable
 Acciaio inossidabile



Jason
CFS-01-195-34 34 cm - 13"

Chest-support with support-rod for adults
 Bruststütze mit Stützstange für Erwachsene
 Soporte de pecho con vara para adultos
 Appui-poitrine avec tige pour adultes
 Supporto per petto con asta per adulti



CFS-01-196-24 24 cm - 9"

Support-rod for children
 Stützstange für Kinder
 Vara para niños
 Tige pour enfants
 Asta per bambini



Lupus
Fig. 1

CFS-01-197-00

for laryngoscopes
für Laryngoskope
para laringoscopios
pour laryngoscopes
per laringoscopi



Fig. 2

CFS-01-198-00

Suction tube for laryngoscopes
Rauchabsaugrohr für Laryngoskope
Tubo de aspiración para laringoscopios
Tube d'aspiration pour laryngoscopes
Tubo d'aspirazione per laringoscopi



Fig. 3

CFS-01-199-00

for laryngoscopes
für Laryngoskope
para laringoscopios
pour laryngoscopes
per laringoscopi



Fig. 4

CFS-01-200-00

Injections drain tube for laryngoscopes
Injektionskanüle für Laryngoskope
Canula de inyección para laringoscopios
Canule à injection pour laryngoscopes
Cannula d'iniezione per laringoscopi